

# Is it an Infectious Disease?

With the **ACON<sup>®</sup>**  
Infectious Disease  
rapid tests, you will  
know in minutes!

- *High sensitivity and specificity*
- *Results in minutes*
- *Room temperature storage*
- *Multiple choices*  
*Format - strip, device or combo device*  
*Packaging - simple or complex*



**KH**

<http://www.kimhung.vn>




**ACON<sup>®</sup>**

Global Diagnostics for Local Markets™

# ACON<sup>®</sup> Infectious Disease Rapid Tests provide accurate results in minutes


## Strip Test

- Immerse the strip into the specimen or diluted specimen
- Read results in minutes




## Device Test


- Add specimen and buffer (if required)
- Read results in minutes




## Whole Blood Strip Test



## Whole Blood Device Test



## Combo Device Test



Product Name	Specimen	Format	Catalog No.	*Read at	Sensitivity	Tests per kit	Dimension
Chlamydia	Swab/Urine	Device	ICH-502✓	10 min.	See Insert	20***	4.0 mm
			ICH-502S✓			25	
Dengue	Whole Blood/Serum/Plasma	Device	IDE-402	10 min.	See Insert	10	4.0 mm
			IDE-402H	10 min.		10	
Gonorrhea	Swab	Strip	IGO-501	10 min.	See Insert	25	4.0 mm
HBV Hepatitis B Virus Combo	Serum/Plasma	Combo Device	IHB-355	15 min.	See Insert	25	4.0 mm
HBsAb Hepatitis B Surface Antibody	Serum/Plasma	Strip	IHBsb-301	15 min.	10 mIU/mL	50	4.0 mm
		Device	IHBsb-302			25**	
HBsAg Hepatitis B Surface Antigen	Serum/Plasma	Strip	IHBsg-301	15 min.	1 ng/mL	50	4.0 mm
		Device	IHBsg-311			25**	
HBsAg Hepatitis B Surface Antigen	Whole Blood/Serum/Plasma	Strip	IHBsg-401	15 min.	1 ng/mL	50	3.3 mm
		Device	IHBsg-411			25**	
HBsAg Ultra Hepatitis B Surface Antigen	Serum/Plasma	Strip	IHBsg-U301	15 min.	0.5 ng/mL	50	4.0 mm
		Device	IHBsg-U311			25**	
HBsAg Ultra Hepatitis B Surface Antigen	Serum/Plasma	Device	IHBsg-U302	15 min.	See Insert	40	4.0 mm
HBsAg Ultra Hepatitis B Surface Antigen	Serum/Plasma	Device	IHBeg-302	15 min.	See Insert	40	4.0 mm
HBsAg Ultra Hepatitis B Surface Antigen	Serum/Plasma	Device	IHBeg-302†	15 min.	See Insert	40	4.0 mm
HBcAb Hepatitis B Core Antibody	Serum/Plasma	Device	IHBcb-302	15 min.	See Insert	40	4.0 mm
HCV Hepatitis C Virus	Serum/Plasma	Strip	IHC-301	10 min.	See Insert	50	4.0 mm
		Device	IHC-302			40	
HCV Hepatitis C Virus	Whole Blood/Serum/Plasma	Strip	IHC-401	10 min.	See Insert	50	4.0 mm
		Device	IHC-402			40	
HIV 1/2 Human Immunodeficiency Virus 1/2	Serum/Plasma	Device	IHI-302	10 min.	See Insert	40	4.0 mm
		Whole Blood/Serum/Plasma	Strip	IHI-401	15 min.	See Insert	
			Device	IHI-411			
			IHI-402			40	
HIV 1/2 Ultra Human Immunodeficiency Virus 1/2	Serum/Plasma	Device	IHI-U302	10 min.	See Insert	40	4.0 mm
HIV 1/2/O Tri-line Human Immunodeficiency Virus 1/2/O	Whole Blood/Serum/Plasma	Device	IHI-T402	10 min.	See Insert	40/10	4.0 mm
			IHI-T402H			10	
H. pylori	Serum/Plasma	Device	IHP-302✓	10 min.	See Insert	40	4.0 mm
		Whole Blood/Serum/Plasma	Strip			IHP-401†✓	
			Device	IHP-402†✓			40
Malaria P.f.	Whole Blood	Device	IMA-402	15 min.	See Insert	40	4.0 mm
Mono Mononucleosis	Whole Blood/Serum/Plasma	Device	IMO-402✓	5 min.	See Insert	20	4.0 mm
Strep A	Throat Swab	Strip	IST-501†✓	5 min.	See Insert	25***	4.0 mm
			IST-501S✓			50	
			IST-501H			1	
		Device	IST-502†✓			20***	
			IST-502F✓			5	
			IST-502S✓			25	
Strep A Twist	Throat Swab	Twist Device	IST-502T†✓	5 min.	See Insert	25***	4.0 mm
Syphilis Ultra	Serum/Plasma	Device	ISY-U302✓	10 min.	See Insert	40	4.0 mm
			Whole Blood/Serum/Plasma			Strip	
		Device				ISY-U411	
			Device			ISY-U402✓	
TB Tuberculosis	Whole Blood/Serum/Plasma	Device	ITB-402	10 min.	See Insert	40	4.0 mm

\*\* Canister packaging, also available in 20 or 30 tests

\*\*\* Includes swabs and/or controls

† Cleared for US 510(k)

✓ CE Marked

## We also offer other rapid diagnostic and medical products for:

Pregnancy (hCG), Ovulation (LH), Menopause (FSH), PSA, FOB, AFP, CEA, Troponin I, Myoglobin, CK-MB, CRP,<sup>NEW!</sup> Drugs of Abuse, Urine Adulteration Strips, Alcohol Saliva Strips,<sup>NEW!</sup> Urinalysis Reagent Strips,<sup>NEW!</sup> Blood Glucose Monitoring Systems, Blood Glucose Visual Strips and others.

**Digital Thermometer,<sup>NEW!</sup> Automatic Blood Pressure Monitor<sup>NEW!</sup> and Urine Analyzer<sup>NEW!</sup> Now Available!**

**To order, please call your local distributor or call ACON for a distributor near you.**

ACON Laboratories, Inc., 4108 Sorrento Valley Boulevard, San Diego, CA 92121, U.S.A. • Tel: 1-858-875-8000 • Fax: 1-858-535-2035 • E-mail: [info@aconlabs.com](mailto:info@aconlabs.com)  
Please visit our website for details: [www.aconlabs.com](http://www.aconlabs.com)