

NEW...

The New Wave of Diagnosis

Dengue Ag + IgG/IgM Rapid Test


*To improve diagnosis during all clinical stages of dengue infection,
Dengue antigens AND antibodies must be tested together!*

Benefits of testing both antigens & antibodies:

- Dengue Antigen is detectable 1 day after the onset of fever and lasts up to 9 days in both primary and secondary infection.
- Primary infection IgM is seen 3-4 days after the onset of infection
- When antibodies are produced, detection of antigen can be inhibited

Duration throughout course of Infection	Ag Test	IgG/IgM Antibody Test	Combination of Dengue Ag & IgG/IgM is the best option for diagnosing Dengue infection during all clinical stages!
1 day after the onset of fever	Begin detectable levels	Negative	
3-4 days after onset of infection	Detectable levels	IgM positive and/or IgG positive during secondary infection	
9 days after the onset of fever	Decline begins	IgG & IgM positive	
Recovery Phase	Negative	IgG positive	

R0062C	Dengue Ag + IgG/IgM Rapid Test	30 tests/kit
R0063C	Dengue Ag Rapid Test	30 tests/kit
R0060C	Dengue IgG/IgM Rapid Test	30 tests/kit
R0061C	Dengue IgG/IgM COMBO Rapid Test	30 tests/kit
E0312	Dengue Ag ELISA kit	96 well



CTK Biotech

6748 Nancy Ridge Dr.
San Diego, CA, USA 92121
858-457-8698 Phone
858-535-1739 Fax
info@ctkbiotech.com
www.ctkbiotech.com